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## Primer Set for *Alicyclobacillus acidoterrestris*

### 【Characteristics】

LAMP(Loop-mediated Isothermal Amplification)method is a novel gene amplification method capturing the following characteristics : ① Only one enzyme is required and the amplification reaction proceeds under isothermal condition<sup>1),2)</sup>, ② It has extremely high specificity because of the use of 4 primers recognizing 6 distinct regions on the target., ③ It has high amplification efficiency and enables amplification within a shorter time., ④ It produces tremendous amount of amplified products which makes simple detection possible<sup>3), 4), 5), 6)</sup>.

This reagent consists of the primer set for the detection of the 16S rDNA *Alicyclobacillus acidoterrestris*<sup>7)</sup>, which is specifically designed for Loopamp DNA Amplification Kit.

### 【Contents of the kit】

48 tests

- |                                      |                  |
|--------------------------------------|------------------|
| (1) Primer Mix. AAT (PM AAT) *1      | 0.12 mL × 1 tube |
| (2) Positive Control AAT (PC AAT) *1 | 0.1 mL × 1 tube  |

\*1 : The notation on each reagent tube is shown in ( ).

### 【How to use】

#### 1. Materials required but not provided

Refer to the package insert for Loopamp DNA Amplification Kit.

#### 2. Preparation of the specimen

Carry out the extraction of DNA from the specimen using a commercially available PrepMan™ Ultra Reagent (Applied Biosystems) or a method equivalent to it.

When handling the sample, adopt appropriate measures to avoid the contamination of *Alicyclobacillus acidoterrestris* to the environment and samples.

#### 3. Reagents preparation

- Take out the reagents stored at  $-20^{\circ}\text{C}$ , and thaw them at room temperature. Once the reagents are thawed, keep them on ice.

- Preparation of master mix. (Operate on ice).

Referring to the table below, dispense the necessary amount of Primer Mix. AAT (PM AAT) of this kit, and 2 × Reaction Mix. (RM), *Bst* DNA Polymerase and Distilled Water (DW) of Loopamp DNA Amplification Kit, into sterilized tube for master mix.

<Reagents>	<Amount : 1 test>	<Amount : 10 tests>
2 × Reaction Mix. (RM)	12.5 $\mu\text{L}$	125 $\mu\text{L}$
Primer Mix. AAT (PM AAT)	2.5 $\mu\text{L}$	25 $\mu\text{L}$
<i>Bst</i> DNA Polymerase	1.0 $\mu\text{L}$	10 $\mu\text{L}$
Distilled Water (DW)	7.0 $\mu\text{L}$	70 $\mu\text{L}$
Total	23.0 $\mu\text{L}$	230 $\mu\text{L}$

#### 4. Operation procedures (Follow the instruction of Loopamp DNA Amplification Kit) .


- Mixing master mix. solution and the sample solution (Operate on ice).
  - Dispense 23  $\mu\text{L}$  of the master mix. into each Loopamp Reaction Tube.
  - Add 2  $\mu\text{L}$  of sample DNA, Positive Control AAT (PC AAT) or Negative Control (DW) into the master mix., and the volume of the solution should be 25  $\mu\text{L}$  in total. Mix the solution well by pipetting or tapping the tube and then spin down. Be careful not to cause air-bubbles when mixing.

#### 2) Amplification and detection

- Incubate at  $65^{\circ}\text{C}$  for 60 minutes.
- After the amplification reaction, use a heat block to inactivate the polymerase and terminate the reaction (5 minutes at  $80^{\circ}\text{C}$  or 2 minutes at  $95^{\circ}\text{C}$ ).

☆ If other final judgment procedures such as fluorescence visual detection or electrophoresis is to be taken, be sure that the enzyme is inactivated. For information on fluorescence visual detection, refer to the package insert for Loopamp Fluorescent Detection Reagent.

☆ For examples of detection cases using this reagent, refer to the products introduction page in Eiken GENOME SITE (URL; <http://loopamp.eiken.co.jp/e/>)

The amplification curve of Positive Control AAT (PC AAT) \*2 (monitored by Loopamp Realtime Turbidimeter)

\*2 : There is no correlation between the initial template number and the turbidity increment time.

### 【Caution for Handling】

- This reagent is designed for research use only.
- If the operator does not have the experience or knowledge in the field of the nucleic acid testing, there's a possibility of false judgment. Therefore, make sure that the reagent is used under the supervision of the experienced and knowledgeable technicians.
- Eiken Chemical Co., Ltd. does not bear any responsibility for false judgment or any consequential damage derived from the false judgment caused by non-capability problems such as operation error.
- Use the reagent before the expiration date, which is labeled on the outer box (Exp.Date).
- If a sample contains more than  $1.0 \times 10^7$  cells/mL of *Alicyclobacillus acidophilus* or *Alicyclobacillus hesperidum*, which are homologous to *Alicyclobacillus acidoterrestris*, it may show cross reaction.
- The reagent tube is made of polypropylene and the main material for kit case is paper. The institution disposing the reagent tube and case should bear the responsibility and abide by the clinical waste disposal regulations, water pollution prevention law, and any other regulation related.

### 【Unit, Storage, Expiration, Code No.】

Product Name	Unit	Storage	Expiration	Code No.
Loopamp® Primer Set for <i>Alicyclobacillus acidoterrestris</i>	48 tests	$-20^{\circ}\text{C}$	1 year	PM0005

### 【References】

- Notomi T. *et al.* : Nucleic Acids Research **28**, No.12, e63 (2000)
- Nagamine K. *et al.* : Clin. Chem. **47**, No.9, 1742-1743 (2001)
- Mori Y. *et al.* : Biochem. Biophys. Res. Commun. **289**, No.1, 150-154 (2001)
- Tomita N. *et al.* : Abstract for the 73rd Annual Meeting of the Japanese Biochemical Society (2000)
- Mori Y. *et al.* : Abstract for the 23rd Annual Meeting of the Molecular Biology Society of Japan (2000)
- Tomita N. *et al.* : Abstract for the 26th Annual Meeting of the Molecular Biology Society of Japan (2003)
- Shimabukuro C. *et al.* : Abstract for the Annual Meeting of the Japan Society for Bioscience, Biotechnology, and Agrochemistry 2004 (2004)