
## < Procedure >

## 【Judgment】

### 1. Judgment of amplification by control reactions

Check that the result of positive control (Control DNA) is [+] and that of negative control (1：1 mixture of Extraction Solution and the embryo washing solution (or the biopsy buffer)) is [-]. If not, conduct the LAMP reaction again as amplification reaction might not proceed properly.

### 2. Judgment of sex

	Reaction Mix. I ( male specific primers )	Reaction Mix. II ( male-female common primers )	Judgment
Result	+	+	Male
	-	+	Female
	+	-	Re-test
	-	-	

## 【Caution for use】

### 1. Storage condition

This reagent kit should be stored at −20℃. To avoid inactivation of reagent components, put the reagent tubes onto ice when reagents are taken out from freezer. Gently spin down reagent if any liquid stays on the wall or cover of its container. Mix thoroughly all reagents prior to use except *Bst* DNA Polymerase. Store reaction tubes in the box provided in order to prevent them from cracking and scratching, which cause false judgment of test result. Store Extraction Solution with its cap tightly closed.

### 2. How to avoid contamination

As the LAMP reaction is extremely sensitive, any contamination of bovine genomic DNA other than the object or amplified product of other specimen causes wrong interpretation of the test result. To avoid such possible contamination, it is strongly recommended to prepare reagents and samples in a specially designated separate room or in a clean bench. Especially, contamination of amplified product of other specimen not only causes false judgment of test result but also pollutes testing environment, and in this case, correct test result may not be obtained unless pollution is completely removed. To avoid this type of contamination, throw away used reaction tube without opening the cap.

### 3. Result

It is recommended to use this test kit under the direction of specialist aboutgenetic testing, as person without such knowledge may misinterpret test result. In any event Eiken Chemical Co., Ltd. shall not be liable for any direct, indirect, special or consequential damages arising out of or in any way connected with the judgment of test results if judgment is done based on the test result that does not originate in the performance of the kit.

### 4. Caution for disposal

The reagent tube is made of polypropylene and the main material for kit case is paper. The institution disposing the reagent tube and case should bear the responsibility and abide by the clinical waste disposal regulations, water pollution prevention law, and any other regulation related.

## 【Unit, Storage, Expiration, Code No.】

Product Name	Unit	Storage	Expiration	Code No.
Loopamp® Bovine Embryo Sexing Kit	24tests	−20℃	1 year	LMP001

## 【References】

- Notomi T. *et al.* : Nucleic Acids Research **28**, No.**12**, e63 (2000)
- Nagamine K. *et al.* : Clin. Chem. **47**, No. **9**, 1742–1743 (2001)
- Mori Y. *et al.* : Biochem. Biophys. Res. Commun. **289**, No.**1**,150–154 (2001)
- Tomita N. *et al.* : Abstract for the 73rd Annual Meeting of the Japanese Biochemical Society (2000)
- Mori Y. *et al.* : Abstract for the 23rd Annual Meeting of the Molecular Biology Society of Japan (2000)
- Tomita N. *et al.* : Abstract for the 26th Annual Meeting of the Molecular Biology Society of Japan (2003)
- Kageyama S. *et al.* : Japanese Journal of Embryo Transfer **25**, No.**3**, 136–140 (2003)
- Hirayama H. *et al.* : Theriogenology **62**, 887–896 (2004)

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
## < サンプルの調製 >

# Bovine Embryo Sexing Kit

This reagent test kit is a LAMP-based sexing of bovine embryo using a portion of fertilized egg. LAMP (Loop-mediated Isothermal Amplification) method is a novel gene amplification having characteristics that ①all reaction can be conducted under isothermal condition using one type of enzyme<sup>1),2)</sup>, ②it has extremely high specificity because of the use of 4 primers recognizing 6 distinct regions on the target, ③it has high amplification efficiency and enables an amplification within a shorter time, and ④it produces extremely large amount of amplified products and is suitable for simple detection<sup>3), 4), 5), 6)</sup>.

#### 【Contents of the kit】

	for 24 tests
(1) Extraction Solution (EX)*	1.5 mL×1 tube
(2) Reaction Mix. I <span> </span> : dNTPs, buffer, male specific primers (RM I)*	0.5 mL×1 tube
(3) Reaction Mix. II <span> </span> : dNTPs, buffer, male-female common primers (RM II)*	0.5 mL×1 tube
(4) <i>Bst</i> DNA Polymerase (Bst DNA Polymerase)*	60 μL×1 tube
(5) Control DNA (Cont DNA)*	0.1 mL×1 tube

\*：The notation on each reagent tube is shown in ( ) .

#### 【Principle】

Amplification of nucleic acids is conducted by the LAMP method using male specific primers recognizing male specific nucleic acid sequence and common primers recognizing male-female common nucleic acid sequence in cells of bovine embryo, and determine sexing by judging whether or not amplification has occurred. Gene amplification is detected using Loopamp End Point Turbidimeter by measuring turbidity of white precipitates of magnesium pyrophosphate, which is a by-product of amplification reaction.

For details of the LAMP method, refer to the Eiken GENOME SITE (URL: http://loopamp.eiken.co.jp/e/).

#### 【Characteristics】

This is a reagent test kit developed in utilizing principle of the LAMP method, and by using Loopamp End Point Turbidimeter specially developed for the LAMP reaction and detection, all reaction from amplification of nucleic acid to its detection can be conducted in a closed system (in the same reaction tube)without necessity of additional step of detection, such as electrophoresis. As presence or absence of amplification is mechanically judged and display of [+] (positive) or [-] (negative) tells sexing of bovine embryo cell, it gives result of sexing easier and in shorter time compared to conventional methods<sup>7), 8)</sup>. The use of male specific primers and male-female common primers for amplification reaction eliminates possibility of misjudgment as much technically as possible. Name of each reagent is indicated on the cap and body of each reagent container in order to avoid mishandling, and the most suitable tube is used for each reagent.

#### 【How to use】

##### 1. Materials required but not provided

- Tools required for micromanipulation
- Loopamp End Point Turbidimeter (Teramecs, Japan)
- Micropipettes (0.5–10 μL, 10–100 μL, 100–1,000 μL)
- Pipette tips with filter
- Aluminum rack for cooling tubes
- Ice (crushed ice) and ice box
- Loopamp Reaction Tube
- Sterilized tube for pretreatment of specimen
- Sterilized tube for preparation of master mix.
- Centrifuge for micro-tubes
- Centrifuge for 8-connected tubes

##### 2. Procedure

###### 1) Pretreatment of specimen

- Prepare necessary quantity of sterilized tubes for pretreatment of specimen and put 6 μL of the biopsied trophectoderm cells of bovine embryo suspended in the buffer into the tubes.
- Add 6 μL of Extraction Solution in each tube to suspend biopsy sample and mix well.
- Leave it stand at room temperature for at least 5 minutes, and spin down after mixing. This refers to Sample Solution.

###### 2) Preparation of master mix. (Operate on ice.)

- Master mix. for male reaction
Pipette necessary volume of Reaction Mix. I (20 μL per specimen)and *Bst* DNA Polymerase (1 μL per specimen) into a sterilized tube and mix well.

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- Master mix. for male-female common reaction
Pipette necessary volume of Reaction Mix. II (20 μL per specimen)and *Bst* DNA Polymerase (1 μL per specimen) into a sterilized tube and mix well.

###### 3) Preparation of reaction mixture and addition of Sample Solution

Use 2 reaction tubes (one for male reaction and the other for male-female common reaction) per specimen. As a single Loopamp Reaction Tube consists of 8 tubes that are connected each other, it can handle four specimens at a time. Dispense 20 μL of Master mix. for male reaction into one reaction tube and 20 μL of Master mix. for male-female common reaction into another tube. Add 5 μL of Sample Solution (total 25 μL) in each tube and spin down after mixing. (Use 5 μL of Control DNA for positive control and 5 μL of 1 : 1 mixture of Extraction Solution and the embryo washing solution (or the biopsy buffer) for negative control instead of Sample Solution.)

###### 4) LAMP reaction

Start the LAMP reaction by setting reaction tubes into the reaction block of Loopamp End Point Turbidimeter and pressing【START/ENTER】button.

###### 5) Detection

Check that 《MEASURE》 sign is shown in the message area on the Loopamp End Point Turbidimeter. Then, transfer reaction tubes from reaction block to detection position and press 【START/ENTER】 button. (Normally, it takes for about 40 minutes from the start of reaction until detection can be done.)

#### 【Caution for operation】

- For biopsy, collect more than 10%(v/v)of trophectoderm. For details of biopsy method and the biopsy washing procedure, refer to the Eiken GENOME SITE (URL: http://loopamp.eiken.co.jp/e/)
- The amount of biopsy sample should be 6 μL. If the volume of the biopsy sample is more than 6 μL, the efficiency of DNA extraction would decrease, which might cause a false judgment.
- Use Loopamp Reaction Tube that is specific for Loopamp End Point Turbidimeter. If other kind of tubes were used as the reaction tube, there would be a possibility for a false judgment due to the optic permeability difference of the measuring light. Check visually if reaction tube has any crack or scratch prior to use. Correct measurement of turbidity cannot be done if a reaction tube has any crack or scratch, and it may cause a false judgment of test result. Take full care when handling reaction tubes.
- Since Extraction Solution and Control DNA are colored, confirm visually if specimen is present in the reaction tube after preparation.
- Extraction Solution gradually deteriorates when exposed to air. Open and close of cap of Extraction Solution should be limited as minimum as possible, and work as quickly as possible, so that the Solution is exposed to air for as short a period time as possible. For the pretreatment of specimen, put the biopsy sample into a tube first, and then add Extraction Solution to it. When storing Extraction Solution, keep its cap tightly closed, and do not aliquot the solution.
- Do not use remaining reagent with other kit even if both kits are the same production lot.
- Handling of the Loopamp End Point Turbidimeter should be conducted according to its operation manual.
- Before starting the LAMP reaction, be sure to check the temperature on the reaction block of Loopamp End Point Turbidimeter is raised to 63° C and that 《READY》sign is indicated on the message area. When the LAMP reaction is over (35 min), the temperature on the reaction block is automatically raised to and maintained at 80° C for 2 min to inactivate DNA polymerase. After the inactivation of the enzyme, the temperature on the heat block is cooled to 63° C.
- When all reactions are over, 《MEASURE》 signal goes on and off in the message area. Take the measurement of the turbidity of the samples while 《MEASURE》 sign is indicated on the message area (for 20 min after all reactions) . Measurement of the turbidity 20 minutes after all reactions may cause false readings, as magnesium pyrophosphate, by-product of amplification reaction, may sediment as precipitates at the bottom of the reaction tube.

#### ■ Protocol for LAMP-based bovine embryo sexing

##### < Preparation of specimen >