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DNA Amplification Kit

This kit is for research use only. Do not use for diagnosis or its assistance for clinical purpose. Read this instruction carefully before use.

【Characteristics】

LAMP (Loop-mediated Isothermal Amplification) method is a gene amplification method capturing the following characteristics : ① The amplification reaction proceeds under isothermal condition^{1), 2)}. ② It has extremely high specificity because of the use of 4 primers recognizing 6 distinct regions on the target. ③ It has high amplification efficiency and enables amplification within a shorter time., and ④ it produces tremendous amount of amplified products which makes simple detection possible^{3), 4)}.

This kit is for amplifying and detecting target gene sequence by the LAMP method in combination with personally designed primers for LAMP, or primer sets sold separately.

【Contents of the kit】

		48 tests	96 tests	192 tests
1.	2 × Reaction Mix. (RM) ^{※1}	0.6 mL ×1	×2	×4
2.	<i>Bst</i> DNA Polymerase (<i>Bst</i> DNA Polymerase) ^{※1}	60 μL ×1	×2	×4
3.	Distilled Water (DW) ^{※1}	1.0 mL ×1	×2	×2
4.	Primer Mix. DNA ^{※2} (PM DNA) ^{※1}	25 μL ×1	—	—
5.	Positive Control DNA ^{※2} (PC DNA) ^{※1}	0.1 mL ×1	—	—

※1 : The notation on each reagent tube is shown in ().

※2 : Primer Mix. DNA and Positive Control DNA are not included in the 96 tests and 192 tests package.

【Principle】

LAMP method is an isothermal nucleic acid amplification method using 4 kinds of primers and DNA polymerase with strand displacement activity. Among 4 kinds of primers, two of them are inner primers, whose 3' end and 5' end are respectively designed to be complementary to two different regions of the target sequence. When the primer-linked strand is replicated, resulting 3' end structure to self-anneals to the complementary region in self-structure, it forms the loop structure at the end.

This stem-loop structure allow the 3' end to initiate self-elongation and another inner primer to anneal to its loop region to synthesize new DNA strand with strand-displacement manner. Through the repetition of these processes, the amplification proceeds, and this enables amplification to continue under isothermal condition with only one kind of enzyme.

For further details of the LAMP method, refer to Eiken GENOME SITE (URL ; <http://loopamp.eiken.co.jp/e/>) .

【How to use】

1. Materials required but not provide

- 1) Sterilized tubes for master mix. preparation (0.5 mL or 1.5 mL)
 - 2) Micropipettes
 - 3) Pipette tips with filter (DNase, RNase free)
 - 4) Loopamp Reaction Tube (sold separately)
 - 5) Aluminum rack for cooling tubes
 - 6) Ice (crushed ice) and ice box
 - 7) Centrifuge for microtube
 - 8) Centrifuge for 8-strip tubes
 - 9) Vortex mixer
- A. For Real-time turbidity detection
Loopamp Realtime Turbidimeter ^{※3}
- B. For Visual fluorescence detection
- 1) Loopamp Fluorescent Detection Reagent (sold separately)
 - 2) Loopamp Realtime Turbidimeter ^{※3}, or incubator with hot bonnet (temperature accuracy within ±0.5℃)
 - 3) Heat block (for enzyme inactivation) ^{※3}
 - 4) UV lamp (wavelength at 240~260nm or 350~370nm) ^{※3}
 - 5) Protective goggle or glass board

※3 : For the information about applicable instrument, reaction termination function and condition for UV irradiation, refer to Eiken GENOME SITE (URL;<http://loopamp.eiken.co.jp/e/>) .

2. Primer design

Appropriate primer design is essential for amplification by LAMP. For the primer design, exclusive primer designing software, PrimerExplorer is available from Eiken GENOME SITE (URL;<http://loopamp.eiken.co.jp/e/>).

When use the highly purified primers, a rapid gene amplification can be performed and stable reproducibility of the amplification can be obtained. The first screening for appropriate LAMP primers might not necessarily require highly purified primers. However, after the primers are determined, it is recommended to use purified FIP and BIP through HPLC or better purification.

3. Reagents preparation

- 1) Take out the reagents stored at -20℃, and thaw them at room temperature. Once the reagents are thawed, keep them on ice. After thawed, tap the reagent tubes gently -- mixing (tapping) -- or invert the tube, or by vortex mixer at about 3 times for 1 second. After mixing well, centrifuge the tubes for a few seconds.

- 2) Preparation of master mix. (**Operate on ice**) .

The following amount of the component is required for one reaction. When used in combination with Loopamp Primer Set, follow the preparation instructions of each primer set.

○ Sample reaction

<Reagents>	<Amount>
Distilled Water (DW)	X μL (ad q.s.)
2×Reaction Mix (RM)	12.5 μL
Primer : FIP	40 pmol
BIP	40 pmol
LF ^{※4}	20 pmol
LB ^{※4}	20 pmol
F3	5 pmol
B3	5 pmol
<i>Bst</i> DNA Polymerase	1.0 μL
Total	23.0 μL/test

※4 : Loop primers are not necessarily required. However, the use of loop primers shortens the amplification time by about one-third⁵⁾.

○ Control reaction

<Reagents>	<Amount>
Distilled Water (DW)	7.0 μL
2×Reaction Mix. (RM)	12.5 μL
Primer Mix. DNA (PM DNA)	2.5 μL
<i>Bst</i> DNA Polymerase	1.0 μL
Total	23.0 μL/test

- 3) Gently tap the tubes for a few times (tapping) , or mix the solution by repeatedly reversing the tube, or mix by vortex mixer 3 times for 1 second. After mixing well, centrifuge the tubes for a few seconds (**Keep on ice**) .

4. Operation procedure

<Mixing of Reagent and Samples> (**Operate on ice**)

Dispense 23.0 μL of master mix. into each Loopamp Reaction Tube.

Add 2.0 μL of sample solution and control into the tubes.

(The volume of the solution should be 25.0 μL in total.)

Negative control : Distilled Water (DW)

Positive control : Positive control DNA (PC DNA)

Mix the solution well by pipetting, or tapping and spin down

(Be careful not to cause air-bubbles when mixing.)

<Amplification>

Set the reaction tubes to the reaction block of Loopamp Realtime Turbidimeter or incubator with hot bonnet, and start the reaction 60~65℃, 30~60 minutes. ^{※5}

Control reaction : 65℃, 60 minutes

Enzyme inactivation (80℃, 5 minutes or 95℃, 2minutes)

(Done by Loopamp Realtime Turbidimeter automatically.)

Turbidity measurement · Judgment

※5 : The reaction condition is dependent upon the characteristics of the primer for use, therefore examine the optimum condition beforehand.

5. Detection

A. For Real-time turbidity detection

Real-time turbidity detection can be carried out with Loopamp Realtime Turbidimeter. For the detailed operation of the equipment, refer to the instruction manual of the equipment.

B. For Visual fluorescence detection

Visual fluorescence detection can be achieved by using Loopamp Fluorescent Detection Reagent⁴⁾. For more information, refer to the package insert for Loopamp Fluorescent Detection Reagent.

6. Amplification plots


Fig1. The amplification plots for control

【Caution for amplification, detection and judgment】

1. LAMP reaction is very sensitive and even the slightest amount of amplified product tainted into the reaction might cause false result. Therefore, avoid this type of contamination and carry out the sample and reagent preparation in different clean benches or separate the area for each operation. In order to prevent contamination, wear gloves and isolation gown, etc. if needed.
2. When this kit used, avoid the contamination of a microbe or nuclease (DNase, RNase).
3. Notice that too much mixing by vortex mixer might inactivate *Bst* DNA polymerase, and assure that vortexing is conducted 3 times for 1 second. The prepared master mix. should be used as soon as possible.
4. Since bubbles in the solution will interfere the turbidity measurement and cause false judgment, try not to cause any bubble when mixing the master mix. and the sample solution. If bubbles are present, spin down to get rid of the bubbles.
5. The cap of the used reaction tubes should not be opened. Contamination of amplified products on other samples may not only cause false judgment of the test result but also pollute testing area. In this case, a correct test result may not be obtained unless pollution is completely removed.
6. Do not carry out electrophoresis of the amplified product.

【Caution for handling (hazard prevention)】

1. The kit is designed for research use only.
2. When handling the sample, always abide by the biohazard counter measures⁶⁾.
3. When ultraviolet lamp is used for the fluorescence visual judgment, do not stare directly at the UV light. Look at it through glass board or protective goggles.

【Caution for handling】

1. This kit should be stored at -20°C . To prevent the reagents from deterioration, only take out the necessary amount of reagents from the freezer before use. (No decline was observed in the kit performance even after repeated freezing and thawing for 20 times in the quality control testing. But, in order to maintain the reagents performance, keep off unnecessary freezing and thawing.)
2. Only use the specified Loopamp Reaction Tube for turbidity or fluorescence detection. Other reaction tubes might have different optical transparency and can cause misjudgment.
3. Take full care when handling reaction tubes, as they are vulnerable to scratches or damages.
4. Check carefully to see if the reaction tubes has any crack or scratch before use. Crack or scratch on the tube might not only cause false judgment but also contaminate the equipment. If the tubes are broken inside the reaction block of Loopamp Realtime Turbidimeter, the reaction solution can spill inside the equipment and cause unrecoverable contamination malfunctioning.
5. Do not expose Loopamp Reaction Tube, master mix. preparation tubes to UV light. A change in color or degeneration caused by ultraviolet lamp sometimes results in misjudgment.
6. If the operator does not have the experience or knowledge in the field of the nucleic acid testing, there's a possibility of false judgment. Therefore, make sure that the kit is used under the supervision of the experienced and knowledgeable technicians.
7. EIKEN Chemical Co., Ltd. dose not bear any responsibility for false judgment or any consequential damage derived from the false judgment caused by non-capability problems such as operation error.
8. Use the kit before the expiration date, which is labeled on the outer box (Exp. Date) .
9. Do not mix different lots.
10. Do not use the container of this kit, and accessories to reuse or other purposes.

【Caution for disposal】

1. Keep the cap of the used tube completely closed and dispose it, according to the relevant regulations and instructions, by incineration or after double bagging it with sealable vinyl bag. To prevent the amplified products from dispersing, **do not conduct autoclave sterilization treatment for disposal.**
2. The reagent tube is made of polypropylene and the main material for kit case is paper.
3. The institution disposing the reagent tube and case should bear the responsibility and abide by the clinical waste disposal regulations, water pollution prevention law, and any other regulation related.

【Storage, Expiration, Unit, Code No.】

Product Name	Storage	Expiration	Unit	Code No.
Loopamp DNA Amplification Kit	-20°C	1 year	48tests	LMP204
			96tests	LMP205
			192tests	LMP206

【References】

- 1) Notomi T. *et al.*, : Nucleic Acids Research **28**, No.12, e63 (2000)
- 2) Nagamine K. *et al.*, : Clin. Chem. **47**, No. 9, 1742-1743 (2001)
- 3) Mori Y. *et al.*, : Biochem. Biophys. Res. Commun. **289**, No.1, 150-154 (2001)
- 4) Tomita N. *et al.*, : Nat Protoc. **3**, No. 5, 877-882
- 5) Nagamine K. *et al.*, : Molecular and Cellular Probes **16**, No.3, 223-229 (2002)
- 6) The guideline for the bio-safety and bio-hazard (by the Japanese Society for Bacteriology): Japanese Journal of Bacteriology **54**, No.3, 667-715 (1999)

Licensed under U.S. Patent #5,814,506.

Manufacturer

EIKEN CHEMICAL CO., LTD.

143 Nogi, Nogi-machi, Shimotsuga-gun, Tochigi 329-0114, Japan

Date of Revision: October 1, 2013 (Ver.5)